Petite&Mars

NEXA

BEZKONTAKTNÍ TEPLOMĚR
BEZKONTAKTNÝ TEPLÓMER
RINTÉSMENTES HŐMÉRŐ
BEZDOTYKOWY TERMOMETR
NON-CONTACT THERMOMETER
TERMOMETRU FĂRĂ CONTACT
BERÜHRUNGSLOSES THERMOMETER
THERMOMÈTRE SANS CONTACT
TERMOMETRO SENZA CONTATTO

NÁVOD K POUŽITÍ
NÁVOD NA POUŽITIE
HASZNÁLATI ÚTMUTATÓ
INSTRUKCJA OBSŁUGI
USER MANUAL
BEDIENUNGSANLEITUNG
INSTRUCȚIUNI DE UTILIZARE
MANUEL D'UTILISATION
MANUALE UTENTE

Výrobce CZ | Výrobca SK | Gyártó HU | Producent PL | Producer EN

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USER MANUAL NON-CONTACT THERMOMETER

Please read this manual carefully before use. If the operating principles are not fully understood, do not use the device. Do not operate the thermometer without first thoroughly reviewing this manual. For further information, contact: ags 92, s.r.o., Drahohejlova 2515/18, 190 00 Praha, Česká republika, info@petiteetmars.com

SAFETY INFORMATION

The device must be used strictly as intended and described in this manual. The manufacturer assumes no liability for damage resulting from improper use. The thermometer is designed to minimize software-related errors through the application of appropriate software engineering standards, risk analysis, and validation processes.

WARNINGS

This thermometer is intended for home use.

Read the instructions, precautions, and technical specifications before use.

This device is designed to measure forehead temperature only. Do not use it for other body parts or purposes.

Do not use a damaged or malfunctioning device.

Allow the device to acclimatize for 30 minutes in environments where ambient temperature has significantly changed (e.g., moving from cold to warm rooms). Operating range: 5°C to 40°C (41°F to 104°F).

Remove the batteries if the device will not be used for an extended period.

The thermometer is not waterproof. Do not immerse in water or other liquids. For cleaning, follow the instructions in the "Cleaning, Care, and Storage" section.

Do not touch the sensor with your fingers.

Avoid measuring temperature in patients with a cold compress on their forehead as it may distort the results.

If ambient conditions change suddenly, ensure no condensation is present on the sensor. If so, clean it according to instructions.

Do not use near strong electrostatic or magnetic fields, which may interfere with accurate measurement.

Do not mix old and new batteries to avoid device malfunction.

Forehead obstructions such as hair, hats, or scarves may affect the result.

Measurement results are indicative. If in doubt, confirm the reading with an alternative method.

PRECAUTIONS

Keep away from children and pets. Store in a dry location, protected from moisture, dust, excessive heat, and direct sunlight. Do not place heavy objects on the device.

Do not dispose of batteries in fire.

Use only the recommended batteries. Do not use rechargeable batteries.

This thermometer is not a substitute for clinical diagnosis.

Do not drop, disassemble, or modify the device.

Do not use if damaged or malfunctioning.

The device contains sensitive components and should be handled with care. Follow storage and operating conditions in the "Technical Specifications."

Do not perform maintenance while the device is in operation.

Do not simultaneously touch the batteries and the patient during use.

Do not use the device if it is damaged or loose, as this may result in injury, incorrect readings, or serious hazards.

No known allergic reactions are associated with use of this device.

Installation and commissioning must be performed in accordance with the accompanying documentation.

INTENDED USE

Petite&Mars non-contact infrared thermometers are intended for periodic and intermittent measurement and monitoring of human body temperature in home and primary care settings.

DEVICE DESCRIPTION**PRINCIPLE OF OPERATION**

The NEXA non-contact thermometer is a portable, battery-operated, reusable device for measuring forehead temperature. It operates using infrared (IR) sensor technology, where the IR sensor emits signals based on object and ambient temperature. An ASIC chip processes these signals into a digital reading, displayed on the LCD screen.

STRUCTURE DESCRIPTION (p. 73) A

- | | |
|----------------------|------------------------------|
| 1. IR Sensor | 4. Power/Scan Button |
| 2. Positioning Light | 5. Battery Compartment Cover |
| 3. LCD Display | |

LCD DESCRIPTION (p. 74) B

- | | |
|--------------------------|---------------------|
| 1. Memory Indicator | 4. Speaker |
| 2. Data Indicator | 5. Battery Status |
| 3. Measurement Indicator | 6. Measurement Unit |

APPLICATION CHART

Measurement Style	Adult	Child
Ear	-	-
Forehead	✓	✓
Non-contact IR	✓	✓

DEVICE SYMBOLS

	Warning		Restriction of Hazardous Substances
	Non-sterile packaging		Batch code

	Refer to operating instructions		Compliance with WEEE Standard
	Operating temperature		DO NOT THROW AWAY – Intended for multiple use
	Operating humidity		Storage and Operating atmosphere pressure
	Manufacturer		Indicates this device is in compliance with MDD 93/42/ EEC, 0598 in the Notified Body Number
	Recyclable		Authorized Representative in the European Community
	Storage temperature		Serial Number
	Storage humidity		IP22 Rating (solid objects ≥ 12.5 mm, water drops at 15°)
	Unique Device Identifier		Medical Device
	Type BF (Body Floating)		

TECHNICAL SPECIFICATIONS

Measurement Unit	°C
Operating mode	Adjusted mode (Body mode)
Measurement site	Forehead
Reference body site	Oral

Rated output range	34,0 °C – 43,0 °C
Laboratory accuracy (Body mode)	34,0 °C – 34,9 °C: $\pm 0,3$ °C 35,0 °C – 42,0 °C: $\pm 0,2$ °C 42,1 °C – 43,0 °C: $\pm 0,3$ °C
Three-color backlight (Color Alarm)	34.0°C–35.7°C: No backlight 35.8°C–37.4°C: Green (normal temperature) 37.5°C–38.5°C: Yellow (slight fever) 38.6°C–43.0°C: Orange (high fever)
Auto Power-Off time	≤ 18 s
Measuring time	Forehead mode: ≤ 2 s Forehead scan mode: 3–10 s
Measuring distance	1–5 cm
Memories	1 reading

Power supply requirements

Batteries	1.5V AAA alkaline x 2 pcs, (IEC type LR03)
Adaptable Range	2.6V–3.6V

Environmental

Operating conditions	Operating temperature: 5°C – 40°C Relative humidity: $\leq 85\%$ RH (no condensing) Atmospheric pressure: 70 – 106 kPa
Transport and storage conditions	Storage temperature: -20°C – 55°C Relative humidity: $\leq 95\%$ RH (no condensing) Atmospheric pressure: 70 – 106 kPa

Dimensions and weight

Weight (without batteries)	47 g
Size	L:122 mm x W:30 mm x H:43 mm

Compliance

Safety Standards	IEC 60601-1 IEC 60601-1-2
------------------	------------------------------

OPERATING INSTRUCTIONS - C

BATTERY INSTALLATION

NOTE: The non-contact infrared thermometer does not operate with depleted batteries and does not support external power supply. Install new batteries before use.

1. Slide the battery compartment cover downward toward the base of the device and remove it.
2. Insert two AAA batteries, ensuring correct polarity (+/-) as indicated inside the compartment.
3. Replace the cover and ensure it clicks into place.

USING THE DEVICE

BEFORE USE: Ensure all warnings and precautions outlined in this manual are read and understood prior to use.

FOREHEAD MODE

1. Press the On/Scan button to power on the device.
2. Aim the thermometer at the center of the forehead (between the eyebrows), maintaining a vertical orientation and a distance of 1–5 cm.
3. Press the On/Scan button again. The positioning light will indicate the measurement spot with a visible dot.
4. The temperature reading will be displayed immediately.
5. Press the On/Scan button to power on the device.
6. Hold the thermometer vertically and close to the forehead.
7. Press and hold the On/Scan button and slowly move the device from one side of the forehead to the other.
8. The positioning light should follow the motion for about 3–10 seconds.
9. The highest recorded temperature during the scan will appear on the screen.

IMPORTANT MEASUREMENT GUIDELINES:

1. If the device has been moved from a colder to a warmer environment, allow it to stabilize in the new setting for at least 30 minutes (5°C to 40°C).
2. The subject's environment should be stable – avoid drafts, fans, or air conditioning directly hitting the individual.
3. If the subject has come from a cooler area, they should remain in the measurement room for at least 5 minutes prior to use.
4. Wait at least 1 second between measurements. After five consecutive readings, allow a pause of 30 seconds before continuing.
5. Do not perform measurements in direct sunlight.

INITIATING MEASUREMENT

1. Turn on the device by pressing On/Scan. The device will perform a self-check and briefly display all segments on the LCD for 1 second.
2. The last recorded temperature will then be shown (Figure 6), followed by a transition to Ready Mode ("----").
3. Forehead Mode: Position the thermometer perpendicular to the forehead at a 1–5 cm distance and press On/Scan. The positioning light will activate and the result will be displayed on the screen.
4. Forehead Scan Mode: Hold the device vertically and close to the forehead. Press and hold On/Scan while scanning slowly across the forehead. After approximately 3–10 seconds, you will hear four beeps ("bip" or "beep-beep"), and the maximum temperature reading will be displayed.

NOTES

After measurement completion and four beeps:

1. Green backlight: Normal temperature (35.8°C–37.4°C)
 2. Yellow backlight: Elevated temperature (37.5°C–38.5°C) – monitor health condition
 3. Orange backlight: High fever ($\geq 38.6^\circ\text{C}$) – cooling and medical consultation recommended
- To maintain accuracy: After five consecutive measurements, wait at least 30 seconds before repeating.

TROUBLESHOOTING

Message	Situation	Solution
Body °C H_i	Temperature taken in not within typical human temperature range (34.0–43.0°C)	The forehead thermometer is for forehead measurement, not other human body site.
Body °C Lo		
Body °C Lo	Measurement over the distance 1- 5 cm	Optimum measurement distance is 1 cm
	Subject's forehead hair, antipyretic sticker head with sweat, etc.	Subject sits quietly with the exposed forehead 5-10minutes before the test.
	Some people's body temperature is lower than the general population.	The main concern is fever temperature.
Body °C Err	Operating temperature exceeds the range of specified temperature.	Move to a room within the operating range and wait 30 minutes before taking temperature.
	The screen flicker, automatic turn off	Replace the battery or the product has been damaged, needs repairs.
	Battery capacity is too low, taking temperature is not possible.	Insert new batteries
POS	Ambient temperature changes too fast	Wait until the ambient temperature will be stable.
	The device does not turn on: 1. Power is off 2. Improper battery installation 3. The battery is exhausted. 4. Display remain blank.	1. Press ON button again 2. Check the battery polarity 3. Replace with a new battery 4. Contact the retailer or service center.

BATTERY REPLACEMENT

1. Open the battery compartment as marked. Ensure the device is off.
2. Remove the old batteries and insert 2 new AAA batteries according to indicated polarity.
3. Close the compartment until it clicks.
4. Do not dispose of used batteries in household waste – return them to a recycling point.
5. If the device does not restart after battery replacement, remove batteries, wait 30 seconds, and reinsert them.

WARNINGS (BATTERY SAFETY)

Do not recharge, dismantle, or incinerate the batteries.

New batteries provide approximately 2000 measurements (average operation time per measurement: 18 seconds).

Use only the recommended battery type.

If the thermometer will not be used for an extended period, remove the batteries.

CLEANING, CARE AND STORAGE

1. The sensor lens is extremely delicate – protect it from contamination and physical damage.
2. Clean only with a clean, soft cloth. Do not use solvents and do not immerse the device in water or any other liquids.
3. Store the thermometer in a clean, dry, dust-free environment.
4. Storage conditions: Temperature range -20°C to 55°C ; relative humidity $\leq 95\%$ (non-condensing).
5. Do not expose the device to direct sunlight, high temperatures, or high humidity, as this may impair functionality.
6. If the thermometer is moved from a cold to a warm environment, allow it to stabilize for 30 minutes before use (operating range: 5°C – 40°C).

DISPOSAL

1. Used batteries must not be disposed of with household waste – dispose of them at designated collection points.
2. At the end of its service life, the device must not be discarded with household waste. Contact local authorities to arrange environmentally responsible disposal in accordance with applicable regulations.

CALIBRATION

1. The device must be installed and used strictly in accordance with the instructions provided in this manual.
2. Warning: Do not use mobile phones or other devices emitting electromagnetic fields near the thermometer – they may interfere with its operation.
3. Warning: The device has been fully tested and verified for proper functionality.
4. Warning: Do not place this device near or stack it with other electronic equipment. If such positioning is necessary, monitor the device to ensure it operates correctly in that configuration.

BODY TEMPERATURE INFORMATION

The human body temperature varies by individual and fluctuates during the day. For this reason, it is recommended that users become familiar with their personal healthy forehead temperature to better assess potential fever conditions.

The typical body temperature range is approximately 35.5°C to 37.8°C.

To determine the presence of fever, compare the current reading with the individual's known normal baseline.

An increase of 1°C or more compared to the baseline generally indicates fever.

Note that measurements from different anatomical sites (rectum, armpit, mouth, forehead, ear) yield different values and should not be directly compared.

Per ASTM standards, the required clinical accuracy of infrared thermometers in the range 37.0°C to 38.9°C is $\pm 0.2^\circ\text{C}$.

WARRANTY

1. The warranty period is 24 months from the date of purchase.
2. The warranty covers manufacturing defects identified within the warranty period.
3. When making a warranty claim, include the product along with a duly completed warranty card and proof of purchase.
4. The product must be used exclusively for its intended purpose.
5. The product must be properly stored and maintained.
6. During the warranty period, the product must be serviced only by authorized service centers – contact the retailer where the product was purchased.
7. The product must be returned clean and in protective packaging.

THE WARRANTY DOES NOT COVER:

1. Defects resulting from failure to comply with the warranty terms.
2. Parts that wear out due to normal use.
3. Damage caused by incorrect connection with other devices.
4. Mechanical damage or damage caused by accidents.
5. Mechanical damage resulting from improper use.
6. Batteries and packaging materials.

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